

Detection of Biotinylated RNA

Either do “dot blot” – pipette biotinylated RNA directly to nylon membrane, followed by UV crosslink as below, or do gel transfer (*in italics below*)

- 1) *Run biotinylated RNA on 0.8 – 1.0% gel in TAE with ethidium bromide.*
 - a. *Run 2 - 5 µg total RNA to detect biotinylation*
 - b. *Run at 200 volts for 15-20 minutes, monitor for heat, can place ice pack under gel box. Larger volume gel box is better for keeping heat low.*
 - c. *Photograph gel to check RNA integrity. Drosophila large rRNA (28S) is heat-sensitive and normally cleaved into two fragments that migrate at same size as small rRNA (18S) – typically see faint large rRNA band and strong small rRNA band*

- 2) *Transfer to Nylon membrane*
 - a. *Use Hybond-N+ membrane (Amersham)*
 - b. *Set up transfer according to Northern protocol of choice*
 - c. *Following transfer, UV cross-link at 1200 joules*

- 3) **Streptavidin-HRP detection**
 - d. Prepare the following solutions:
 - i. Blocking Soln. 500 ml
 - NaCl: 3.65 g
 - Na₂HPO₄: 1.20 g
 - NaH₂PO₄: 0.5 g
 - SDS: 50.0 g

*need to heat soln. to 37° C for SDS to dissolve
 - ii. Wash I: 1:10 dilution of Blocking soln.
 - iii. Wash II: 500 ml of 10X stock
 - Tris: 6.0 g
 - NaCl: 2.9 g
 - MgCl₂: 1.0 g

pH to 9.5 w/ HCl, use at 1:10 dilution
 - e. Block membrane in Blocking Soln,
 - i. 30 min.
 - f. Incubate with Streptavidin-HRP (Pierce) diluted 1:5000 in Blocking Soln
 - i. 5 min.
 - g. Wash in Wash I:
 - i. 20 min., 2X
 - h. Wash in Wash II:
 - i. 5 min., 2X
 - i. Add ECL reagent for 1 min. and expose to film / measure chemiluminescence on Chemidoc